**Bedienungsanleitung  
Instruction Manual  
Mode d'emploi  
Instrucciones de uso**

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Kaufdatum:.....

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## CAUTION!

To work with this microscope, sharp and pointed aids are being used. Please take care that this microscope and its accessories are stored at a place out of reach of children. Let children only work with this microscope under an adult's supervision! Keep packing material (plastic bags etc.) away from children!

### 1. Package contents

(several accessory parts are premounted)

Microscope head  
10x WF eyepiece  
Objective (4x, 10x, 40x, 100x oil)  
Condenser  
Filter plate (blue)  
AC adapter  
Dust proof cover

### 2. Components

- 1 Eyepiece
- 2 Eyepiece holder
- 3 Microscope head
- 4 Objective revolver
- 5 Objective
- 6 Slide clamp
- 7 Condenser
- 8 Iris diaphragm (with adjusting lever)
- 9 Filter holder
- 10 Condenser height adjustment knob
- 11 Illumination
- 12 Coarse-adjusting wheel
- 13 Fine-adjusting wheel
- 14 Lengthwise knob of mechanical desk
- 15 Cross knob of mechanical desk
- 16 Dimmer
- 17 On/Off switch
- 18 Status-LED (red, green)
- 19 AC power socket

### 3. General/Location

Before you set up the microscope, you must choose a suitable location.

First you must make sure that your microscope is on a stable and solid surface.

### 4. Electric illumination

It works, by plugging the connection cable into the microscope socket (19) and plugging the AC adapter into the mains socket. The switch at the rear of the microscope must be switched on.

HINT for the rechargeable battery:

Your microscope is fitted with one built-in rechargeable battery. A fully charged battery allows ca. 48 hours of observation. Before using the microscope for the first time the rechargeable battery should be charged for 16 hours. Please charge it by connecting the microscope socket (19) with the included AC adapter and plugging the AC adapter into mains socket.

### 5. Observation

You begin each observation with the lowest magnification (eyepiece 10x and objective 4x). So, the object to be viewed can easily be centred and properly focussed. The higher the

magnification the more light is required for good image quality. For applications with the 100x lens, it is recommended that examinations are done only with oil immersion.

Begin with a simple observation. Turn the objective revolver (4) onto the smallest magnification and insert the 10x eyepiece. Now place a preparation directly under the objective on the microscope table. Move the preparation on the table with the two adjusting screws (14 and 15).

Now look through the eyepiece (1) and turn the coarse adjustment (12) carefully until the image becomes visible. The optimal sharpness is achieved by the fine adjustment (13).

**Make sure that you never overtighten the fine adjustment.**

### 6. Condenser

Use the condenser height adjustment knob (10) to properly bundle the light for your object.

### 7. Iris diaphragm

An iris diaphragm (8) is available if you want better image quality. This makes it possible to change perspective depth. The smaller the aperture the greater the perspective depth. However, image brightness is also reduced. You can compensate for this by adjusting the brightness using the regulator (16).

### 8. Filter holder

Push the filter holder (9) forwards and then insert the blue filter if you wish to use it. The filter holder can then be pushed back again. This filter is used to show natural colours.

### 9. Oil immersion

The immersion oil is only necessary for the 100x objective. Give one drop of immersion oil onto the upper condenser lens. With the slide clamp (6), fix your prepared slide on the mechanical desk. Using the condenser height adjustment knob (10), move up the condenser (7) carefully, until the slide touches the immersion oil. Focus the preparation with the 100x objective and the fine focuser as good as possible. Now turn the objective revolver a little bit to the side and give one drop of immersion oil onto the preparation's covering glass.

Turn back the objective revolver. The objective must dive into the oil! Now the preparation can be focused. Don't forget to clean objective and preparation with alcohol after observing.

### 10. Care and cleaning

1. Avoid disassembling the microscope!

Because of being a precise optical instrument, do not disassemble the microscope casually. That may cause serious damage to its performance. Do not disassemble the objectives.

2. Avoid contamination!

You should make sure that dust or moisture do not come in contact with your microscope. Avoid putting fingerprints on any optical surfaces! If dust or dirt nevertheless have gotten onto your microscope, or the accessories, remove it as described below.

3. Cleaning non-optical parts

Disconnect the plug from mains socket before cleaning

non-optical parts of the microscope (e. g. the frame or the stage). Firstly, remove dust or dirt from your microscope, or accessories, using a soft brush. Then clean the soiled area with a soft, lint-free cloth; the cloth may be moistened with a mild detergent solution. Be sure that the microscope is dry before using. For cleaning the painted parts or plastic components, do NOT use organic solvents such as alcohol, ether, acetone, xylene or other thinner etc.!

4. Cleaning optical parts

In order to enhance the optical quality, the lenses of eyepieces or objectives are coated. They should not be wiped because dry dirt or dust may scratch the coating. It is best to remove the parts to be cleaned from the frame prior to cleaning. Always blow loose dust away first. Use lens tissue of good quality or soft fabric moistened with a small amount of pure alcohol (available in the chemist's shop); wipe the lens surface clean.

5. Cleaning the oil immersion lens

After working with oil immersion, the immersion oil should be removed from the lens of the 100x objective at the end of each workday. In order to do this, use a lens tissue moistened with a small amount of pure alcohol (available in the chemist's shop). Don't forget to clean the preparation, too.

6. Being not used

After usage, cover the microscope with the dust cap (made of plastic) and place it where it is dry and mouldless. In the case of non-use for a longer time, you should stow the microscope and the accessories in their correct containers again. We suggest the storage of all objectives and eyepieces in a closed container with drying agent.

**Remember:**

**A well maintained microscope will keep its optical quality for years and thus maintain its value.**

### 11. Technical data

Monocular head  
Mechanical desk (with vernier alignment)  
Eyepiece: DIN WF 10x  
Objectives: DIN 4x / 10x / 40x / 100x oil  
Magnification: 40x -1000x  
Illumination: Accu-LED  
Accessories: Immersion oil

**Technical alterations and errors excepted.**

**Aperture:**

The numerical Aperture (the unit of measurement for the resolution factor of the objectives) is engraved beside the magnification-sign on the objectives.

**The magnification for microscopes**

is calculated by formula:  
Eyepiece magnification x objective magnification factor = total magnification  
Example: 10x eyepiece x 100x lens = 1000x total magnification

### 12. EEC Conformity explanation

Meade Instruments Europe GmbH & Co. KG, resident in 46414 Rhede/Westf., Gutenbergstr. 2, Germany, explains the agreement with in the following specified EEC guidelines for this product:

EN 61558-2-6: 1997  
EN 61558-1: 1997 +A1

Product description: Biological-type microscope  
Model: BRESSER Erudit DLX

Rhede, 20th July 2007

Meade Instruments Europe GmbH & Co. KG

Helmut Ebbert  
Managing director

### 13. Warranty

The period of warranty is 2 years, beginning on the day of purchase. Please keep the cash receipt as evidence of purchase. Devices which become defective during the warranty period can be returned to the dealer where the device was bought. The repaired device or a new one will then be returned to you. In the case of defects which occur after the end of the warranty period, the devices can also be returned. However, repairs which become necessary after the end of the warranty period will be subject to a service fee.

**Important:**

**Make sure to return the device carefully packed in the original packaging in order to prevent transport damage. Please also enclose the cash receipt (or a copy). This warranty does not imply any restriction of your statutory rights.**

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